

Laboratory 7

Measurement of magnetic susceptibility and Barkhausen Effect

<https://vlab.amrita.edu/?sub=1&brch=192&sim=854&cnt=1>

<https://vlab.amrita.edu/?sub=1&brch=192&sim=347&cnt=1>




Ajay Nath

Quincke's Method

Aim: - To determine the volume magnetic susceptibility of Manganese sulphate solution at different concentrations.

- Used to determine the magnetic susceptibility of a paramagnetic substance in the form of liquids, aqueous solutions and liquified gases.
- This method is based on the force experienced by a magnetized material in a non-uniform magnetic field.
- The solution under investigation is placed in a vertical U tube which is then placed between electromagnets.
- When current is switched ON, strong field is experienced by upper surface of the column and weak field is experienced by lower surface.



The magnetic susceptibility χ_m is a proportionality constant which is dimensionless and indicates the degree of magnetization of a material in response to an applied magnetic field.


The susceptibility of the paramagnetic solution can be given by:

$$\chi_{sol} = \frac{2gh(\rho - \sigma)}{\mu_0 H_m^2}$$

Setup 1

Step 1: Insert the probe.

Quincke's Method



VARIABLES

Select set up:

Magnetic field Vs Current

REMOVE PROBE

Current: 2 A

RESET


RESULT

Step 2:
Adjust the
current.

Step 3:
Note the
magnetic
field.

Setup 2

Quincke's Method



Step 1: Adjust the focus and microscope position.

Step 2: Set the molarity and vary the current.

Step 3: Note the magnetic susceptibility of the solution.

VARIABLES

Select set up:

Quincke's setup

Adjust Focus Knob:

Current: 2 A

Molarity: 0.5 M

Adjust Microscope:

☒ Show Result


RESET

RESULT

Magnetic Susceptibility of the solution:
 9.13×10^{-6}

4. Now Select **Quincke's setup**

Press **Esc** to exit full screen



VARIABLES

Select set up:

Quincke's setup ▼

Adjust Focus Knob:



Current: 1 A



Molarity: 0.5 M



Adjust Microscope:



☒ Show Result

RESET

RESULT

Magnetic Susceptibility of the solution:
 9.13×10^{-6}

Further steps to follow: -

ng the

current
change

the
level

of
adjust
the
solution

Now
cope
at
different

focus
knob
values

adjust
the
mola

Plot a
note
graph
between

Magn
netic
molar

ity
and
Magn
etic
Susce

ptibili
ty of
the
soluti

Observation Tables

Setup 1

| Current I (A) | Magnetic Field B (Gauss) |
|---------------|--------------------------|
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |

**Plot a graph between current and magnetic field

Setup 2

Current = A


| Molarity (M) | Magnetic Susceptibility (χ_{sol}) $\times 10^{-6}$ |
|--------------|--------------------------------------------------------------|
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |

**Plot a graph between molarity and susceptibility.

Barkhausen Effect:


https://en.wikipedia.org/wiki/Barkhausen_effect

- Series of sudden changes in the size and orientation of ferromagnetic domains, or microscopic clusters of aligned atomic magnets, that occurs during a continuous process of magnetization or demagnetization.



Procedure: -

1. We can firstly remove the rod and look at the output graph by moving the magnet.
2. Then we can insert the rod of ferromagnetic material and look at output graph of EMF produced with time.
3. We can change the no. of coils and see how produced EMF changes.



Domains in ferromagnetic material

Practical use: -

The amount of Barkhausen noise for a given material is linked with the amount of impurities, crystal dislocations, etc. and can be a good indication of mechanical properties of such a material.

**Thank
You**